



Electrodialysis Cell Unit

PCCell ED 64 0 02



Technical Data



PCCell GmbH
Lebacher Str. 60
66822 Heusweiler
Germany
Phone: ++49 - 6806 - 603732
Fax: ++49 - 6806 - 603731
E-mail: pccell@electrodialysis.de
www.pccell.de

Version: Jan 06

Table of Contents

| | |
|--|----|
| 1. General..... | 2 |
| 2. Function..... | 3 |
| 3. Cell and Parts..... | 5 |
| 4. Technical Data..... | 8 |
| 5. Application Examples..... | 10 |
| 6. Maintenance and Replacement..... | 11 |
| 7. Application areas..... | 12 |
| 8. General Remarks / Safety..... | 13 |
| 9. Warranty | 13 |
| 10. Further Information / Contact Address..... | 14 |

1. General

The electro dialysis cell unit PCCell ED 64 is used in laboratory electro dialysis processes to remove ions from one solution (diluate). The ions are collected in another solution (concentrate). The electro dialysis unit PCCell ED 64-2 allows to carry out different types of experiments for a variety of applications, to examine the characteristics of ion exchange membranes in use. It is conceived as an easy-to-manage laboratory cell.

2. Function

The electro dialysis unit *PCCell ED64-2* consists of an anode chamber, a cathode chamber and a membrane stack between them. With this constitution, a variety of experiments, like desalination, can be processed.




Fig. 1: Functional setup of an ED stack. Salts are removed in cells called "Diluate" and are collected in the Concentrate. Beside this, the electrodes need a solution, the Catholyte and the Anolyte.

To run a standard ED, the membrane stack consists of n (typically 5, 10, 50 or even 100) cell pairs, which are formed by $n+1$ cation exchange membranes, n anion exchange membranes and $2n$ spacers.

At the shown polarity (Fig. 1), one of the cell systems is the diluate and the other one is the concentrate. If the polarity is changed, the function of the cell system changes accordingly.

A complete ED System is set up by this ED Cell in combination with an ED pump unit, e.g. PCCell B-ED 1 and the external solvent tanks build (Fig. 2).




Fig. 2: A complete ED setup: it consists of the stack (upper dotted rectangle), the ED pump unit (middle) and the external electrolyte containers (below).

Please remark in this context, that the ED cell is only one part of the complete system and will work properly only in combination with the other parts.

3. Cell and Parts

The PCCell ED 64 consists of several parts. They are described as follows:




Fig. 3: PCCell ED 64 setup and functional parts

| # | Description | Art No. |
|---------|----------------------|---------|
| 1,2 | Electrode end plates | |
| 3, 4, 5 | Membrane stack | |
| 4 | Screw set | |




Fig. 4: PCCell ED 64 0 02 Electrode End Plates.

| # | Description |
|----|--|
| 1 | Electrode Block (PP) |
| 2 | Electrolyte liquid inlet |
| 3 | Concentrate (Diluate) inlet |
| 4 | Clamping Bolt hole |
| 5 | Concentrate (Diluate) connector hole to spacer |
| 6 | Electrode |
| 7 | Membrane and Spacer size and position |
| 8 | Diluate (Concentrate) connector hole to outlet |
| 9 | Electrode compartment |
| 10 | Diluate (Concentrate) outlet |
| 11 | Electrolyte liquid outlet |




Fig. 5: PCCell ED 64 0 02 Membrane size. Size given is in mm. Dotted square area in the middle of the membrane is the active membrane area of the cell.

4. Technical Data

Stack

| | |
|----------------------|--|
| Membrane size | 110 x 110 mm |
| Active membrane area | 64 cm ² |
| Membrane spacing | electrode - membrane: ca. 1 mm over cells: 0,5 mm |
| Number of membranes | Maximum 25 cell pairs |
| Current Connectors | 4 mm banana plugs |
| Processing length | 80 mm |

Medium Contacting Materials End plate materials

| | |
|------------|-----------------------|
| Cellframe | polypropylene |
| Tubes | polyethylene |
| Electrodes | Titane, Pt/Ir coating |

Spacer options

| | PVC | Silicone | Viton | EPDM | Food approved silicone / polyethylene | Silicone / Polyester | Silicone / PVDF | Silicone / Polyethylene |
|-------------------|-----|----------|-------|------|---------------------------------------|----------------------|-----------------|-------------------------|
| End spacer | 1 | 1 | 1 | 1 | | 1 | 1 | |
| Spacer ED 64 0 04 | | | | | 1 | 1 | 1 | 1 |
| Spacer ED 64 0 02 | | | | | 1 | 1 | 1 | 1 |

1 available

Dimensions and Weight (ca)

| | |
|--------|--------|
| width | 165 mm |
| depth | 150 mm |
| height | 190 mm |
| weight | 2.5 kg |

Electrical Connecting data

| | ED 64 |
|--------------|---|
| Type | only connect the cell with a galvanically isolated DC current circuit. DANGER! |
| Max Amp. | 5 (depend on application, temperature etc.) |
| Max. Voltage | Max. 2 V / cell pair |

Hydraulical Connecting data

| | ED 64 |
|--|---|
| Type | connectors for 8 mm id tubes |
| Flow through electrode circuits: | nominal 150 l / h |
| Nominal flow through concentrate and diluate per single cell | 4-8 l / h (10 cell pairs result in 40-80 l / h Flux through concentrate or diluate) |
| Max. pressure | transmembrane pressure has to be kept zero: Never pump only one of the diluate / concentrate circuits alone! |
| pressure drop over cell | Max. 0,5 bar |

5. Application Examples

The PCCell ED 64 is intended to be driven in a batch process. Process length and nominal flowthrough are given above.

Figure 7 shows an example of a batch desalination (conductivity of diluate against time). The effect of a single pass desalination at the start and stop time result in a conductivity jump. It depends on the current, flowing and other factors. The plot shows, that it is - more or less - proportional to the current.




Fig. 6: A batch desalination of aqueous NaCl, about 14 l / m².

Figure 6 shows the salt (calculated as NaCl) removal in dependence of the current at theoretical current efficiency (c_e) and at 85% c_e . With the PCCell ED 200 you can expect for sodium chloride c_e 's in the range between 90 and 95 %. It depends on current density, concentration and other factors. The amount is given per cell pair. A 25 cell pair- unit will make 25 times of this.




Fig. 7: Transport rate of one ED cell pair in dependence of the Amperage.

6. Application areas

Desalination of salt water

Stabilisation of wine

Whey demineralisation

Pharmaceutical application

Pickling bath recycling

7. General Remarks / Safety

By running an electro dialysis with this cell unit, concentrated acids and bases, which are corrosive, may be produced. Adequate protective measures have to be taken. At the electrodes, explosive gases and aerosol may be produced. Also in this case, appropriate protection has to be ensured. The cell has to be run in a tank large enough to collect any liquids passing out of the system.

8. Further Information / Contact Address

For further information, visit our web site www.electrodialysis.info. In case of any technical questions, please contact

| |
|---|
| <p>PCCell GmbH Dr. Patrick Altmeier Phone: ++49-(0)6806/603730 Fax: ++49-(0)6806/603731 email: pccell@electrodialysis.de</p> |
|---|

Note:

The information in this handling instructions is presented in good faith, and all recommendations or suggestions are made without guarantee. The products are intended for use by persons having technical skill, at their own discretion and risk. PCCell is not responsible for any risks or liabilities which may result from the use of its products.