

Cellulose Nitrate (CN) Disc Membrane:

Major Advantages:

- BSA protein binding capacity of 160 ug/cm²
- High wick rates (lateral flow)
- Uniform pore structure
- Superior strength
- Exceptional sensitivity
- Low extractables
- Lot-to-lot consistency

Ordering Information

Pore size (µm)	Diameter (mm)					
	13	25	47	90	142	293
	100 per package			25 per package		
0.1	CN01013	CN01025	CN01047	CN01090	CN010142	CN010293
0.22	CN02213	CN02225	CN02247	CN02290	CN022142	CN022293
0.45	CN04513	CN04525	CN04547	CN04590	CN045142	CN045293
0.65	CN06513	CN06525	CN06547	CN06590	CN065142	CN065293
0.8	CN08013	CN08025	CN08047	CN08090	CN080142	CN080293