

## Nitrocellulose membrane for lateral flow



### **Description:**

Nitrocellulose membrane is specially developed for use in rapid tests. It is a pure nitrocellulose membrane supported on PET film and is suitable for all lateral flow applications.

### **Features:**

- 100% pure nitrocellulose
- 5 µm pore size
- High binding capacity
- Fast speed
- Low background
- High sensitivity
- Supported on PET film
- Good for most of rapid tests – pregnancy, drugs of abuse, infectious diseases

**Other sizes are available upon request.**