

Syringe Filter Holders



For filtering small volumes of fluid, especially where aseptic conditions must be maintained, one of the most useful devices is a membrane filter holder that attaches to a syringe. A variety of different syringe filter holders are manufactured by Membrane Solutions.

Both re-usable and disposable units are available, and the disposable ones are marketed pre-sterilized with filter in place. For ultracleaning applications, for filtering small volumes of liquids, and for processing clinical specimens for viable counting, syringe adapters are extremely useful. Any laboratory doing aseptic work should have a supply of these units on hand for routine or emergency use.